

Culturing Izon ES Cell Clones

Cell Line Information

Izon cell clones are feeder-dependent.

Parental ES cells: * W9.5, isolated from 129S1/Sv mouse strain
 * Bruce4, isolated from C57BL/6 mouse strain
 Injection of cells into C57BL/6 blastocysts will produce agouti chimeras. These cell lines are mycoplasma free.

Feeder cells: MEF feeder cells express very small amounts of Leukemia Inhibitory Factor (LIF) and neomycin phosphotransferase (*Neo*). This cell line is mycoplasma free.

Reagents and Supplies

<u>Item</u>	<u>Vendor</u>	<u>Catalog Number</u>
DMEM, high glucose, with glutamine	Gibco	11965-092
Penicillin/Streptomycin, 100 u/ml	Gibco	15140-122
L-Glutamine	Gibco	25030-081
Sodium Pyruvate	Gibco	11360-070
NE Amino Acids	Gibco	11140-050
Trypsin EDTA (1X, 0.25%) (supplement the 0.25% with 1% Chick Serum)	Gibco	25300-054 (0.05%) 25200-072 (0.25%)
LIF (ESGRO)	Gibco	13275-029
Fetal Bovine Serum, Defined	Hyclone	SH30070.03
PBS (1X without Ca or Mg)	Gibco	14190-144
DMSO, 100 ml (2X Freezing medium: FBS with 20% DMSO - make fresh as required)	Sigma	D2650
Mitocycin C 10x2mg (10µg/ml Inactivation media: add 2mg of Mitomycin C to 200 ml MEF Feeder medium -may be stored at 20°C for up to 6 months)	Sigma	M0503
2(β)-Mercaptoethanol (1000x working soln: add 70µl 2-Mercaptoethanol to 9.93ml PBS. Store at 4°C, and make fresh every 2 weeks)	Sigma	M-7522
Gelatin, 2% (0.1% working soln: add 25 ml of 2% solution to 475ml of PBS. Store at 4°C)	Sigma	G1393

MEF Feeder Medium (sterile filter through 0.2 μ M filter unit)

<u>Reagent</u>	<u>Stock Conc.</u>	<u>Final conc.</u>	<u>Total: 500ml</u>
DMEM		1x	435ml
FBS	100%	10%	50ml
L-Glutamine	200mM	2mM	5ml
Sodium Pyruvate	100mM	1mM	5ml
Pen/Strep	10,000U/ml	100U/ml	5ml

Izon / Bruce4 Cell Medium (sterile filter through 0.2 μ M filter unit)

<u>Reagent</u>	<u>Stock Conc.</u>	<u>Final conc.</u>	<u>Total: 500ml</u>
DMEM		1x	403.5ml
FBS	100%	15%	75ml
L-Glutamine	200mM	2mM	5ml
Sodium Pyruvate	100mM	1mM	5ml
NE Amino Acids	100mM	1mM	5ml
Pen/Strep	10,000U/ml	100U/ml	5ml
LIF	10 ⁷ U/ml	1000U/ml	1.0ml
1000x β ME	5.5 x 10 ⁻⁶ M	1 μ M	0.5ml

MEF Inactivation Media (sterile filter through 0.2 μ M filter unit)

<u>Reagent</u>	<u>Stock Conc.</u>	<u>Final Conc.</u>	<u>Volume</u>
MEF Feeder medium	1x	1x	200ml
Mitomycin C	2mg powder	10ug/ml	2mg
Total Volume			200 ml

Preparing MEF Feeder Cells

Thawing, Expanding and Treating Active MEF Cells

1. Thaw 1 vial of MEF cells (approximately 1.5-2 x 10⁶ cells/vial) in a 37°C water bath and dilute into 10 ml of pre-warmed MEF feeder medium.
2. Pellet the cells by spinning for 4 minutes at 1000 rpm in a bench-top clinical centrifuge.
3. Aspirate off medium and gently resuspend cells in 5 ml of pre-warmed MEF feeder medium.
4. Transfer cell suspension to a 6 cm gelatinized dish, and grow at 37°C in a humidified 5% CO₂ incubator (should be confluent within 2-3 days).
5. When confluent, aspirate medium off and wash with 5 ml of pre-warmed PBS, pipetting it away from the cells. Rock dish gently and aspirate medium. Repeat.
6. Cover cells with 1 ml of 0.05% trypsin solution and incubate at 37°C for 4 minutes or until cells are uniformly dispersed into small clumps.
7. Add 5 ml of MEF Feeder medium to inactivate the trypsin, and pipette vigorously to make single cell suspension (we recommend 15 times).

8. Spin for 4 minutes at 1000 rpm.
9. Aspirate off medium and gently resuspend cells in 20 ml of pre-warmed MEF Feeder medium.
10. Split the cell suspension onto two gelatinized 10 cm tissue culture dishes, and grow at 37°C in a humidified 5% CO₂ incubator (should be confluent within 2-3 days).
11. To mitotically inactivate, replace medium with 10 ml Inactivation medium (2mg of Mitomycin C to 200 ml MEF feeder medium), and incubate in a 37°C humidified 5% CO₂ incubator for 2.5 hours. Aspirate Inactivation medium, and rinse three times with pre-warmed PBS; aspirating completely between rinses. These dishes are now ready to use.
12. If you wish to freeze the cells for later usage, trypsinize and pellet the cells as before, but with 1.5 ml of 0.05% trypsin solution, and inactivate the trypsin with 5 ml medium.
13. For each 10 cm dish, count cells, and resuspend in an equal volume of MEF Feeder medium and 2X Freezing medium; to a density of 1.5-2 x 10⁶ cells/0.5 ml. Decant 0.5 ml aliquots into labeled cryovials.
14. Immediately place cryovials in a Styrofoam container or temperature controlled freezing vessel.
15. Freeze vials in a -80°C freezer. After 24 hours, transfer cryovials to liquid or vapor-phase nitrogen for longer term storage.

Plating Mitotically Inactive MEF Feeder Cells

1. Coat a 6 cm tissue culture dish with 0.1% gelatin and aspirate off immediately before use.
2. Thaw 1 vial of mitotically inactive MEF Feeder cells (approx. 1.5-2 x 10⁶ cells) in a 37°C water bath and dilute into 10 ml of pre-warmed MEF feeder cell medium.
3. Pellet the cells by spinning for 4 minutes at 1000 rpm.
4. Aspirate off medium and gently resuspend cells in 5 ml of pre-warmed MEF Feeder cell medium.
5. Transfer cell suspension to the 6 cm gelatinized dish, and grow at 37°C in a humidified 5% CO₂ incubator.
6. The feeders are ready for use after 6-12 hours, or may be maintained in the incubator for a maximum of 8-10 days.

Thawing Izon ES Cell Clones

1. Thaw 1 vial of ES cells (approximately 5 x 10⁶ cells/vial) in a 37°C water bath and dilute (drop wise) into 10 ml of pre-warmed Izon/Bruce4 cell medium.
2. Pellet the cells by spinning for 4 minutes at 1000 rpm.
3. Aspirate off medium and gently resuspend cells in 5 ml of pre-warmed Izon/Bruce4 cell medium.

4. Aspirate the old medium from your 6 cm mitotically inactive MEF Feeder dish.
5. Transfer the ES cell suspension to the feeder dish, and grow in a 37°C humidified 5% CO₂ incubator (Important: 10% CO₂ incubator for Bruce4-derived cell lines).
6. Change medium the following day to remove dead cells and residual DMSO.
7. Change medium daily until 80% confluent (approx. 1.5-2 x 10⁷ cells); should take 2-3 days.
8. When confluent, the 6 cm dish may be split in two; half for microinjection and half to expand for freezing.

Expansion of Izon ES Cell Clones for Microinjection and Future Use

1. 1 day prior, prepare one 10 cm MEF Treated Feeder dish, using the plating density guide below.
2. The next day, aspirate off the old medium prior to plating ES cells.
3. On the day, wash the confluent 6 cm ES cell dish once with 5 ml PBS.
4. Cover the cells with 1 ml of 0.25% trypsin solution and incubate at 37°C for 4-5 minutes or until cells are uniformly dispersed into small clumps.
5. Add 5 ml of Izon/Bruce4 cell medium; to inactivate the trypsin, and pipette vigorously to make single cell suspension (we recommend 10-15 times).
6. Split the cell suspension in half, placing 2.5 ml each into 15 ml centrifuge tubes (labeled 'Expansion' and 'Microinjection').
7. Spin both tubes for 4 minutes at 1000 rpm.
8. For the 'Expansion' cells; aspirate off the supernatant and resuspend the pellet in 10 ml Izon/Bruce4 cell medium. Transfer the cell suspension onto the 10 cm MEF Treated Feeder dish prepared the day before. Grow in a 37°C humidified 5% CO₂ incubator (Important: 10% CO₂ incubator for Bruce4-derived cell lines). Change medium daily until 80% confluent (should take 2-3 days).
9. For the 'Microinjection' cells; aspirate off the supernatant and resuspend the pellet in 150-400 µl microinjection medium (Hepes Buffered DMEM with 5% FBS; filtered through 0.2 µm filter unit, can be aliquoted and stored at -20°C for up to 12 months). Immediately place the cell suspension on ice, and microinject within 1-2.5 hours.

Freezing Expanded Izon ES Cell Clones

1. Wash the confluent 10 cm ES cell dish once with 10 ml PBS each.
2. Cover the cells with 1.5 ml of 0.25% trypsin solution and incubate at 37°C for 4-5 minutes or until cells are uniformly dispersed into small clumps.

3. Add 5 ml Izon/Bruce4 cell medium to inactivate the trypsin, and pipette vigorously to make single cell suspension (we recommend 10-15 times).
4. Spin for 4 minutes at 1000 rpm.
5. Aspirate supernatant and resuspend the pellet in an equal volume of Izon/Bruce4 cell medium and 2X Freezing medium (we would recommend 8-10 vials containing 0.5 ml aliquots; per 10 cm dish). Decant into labeled cryovials.
6. Immediately place cryovials in a Styrofoam container or temperature controlled freezing vessel.
7. Freeze vials in a -80°C freezer. After 24 hours, transfer cryovials to liquid or vapor-phase nitrogen for longer term storage.

MEF Treated Feeder Cells – Plating Density Guide

Plate/Dish	Feeder Cell Density
96 well	1.5×10^4
24 well	1×10^5
3.5 cm dish	8×10^5
6 cm dish	$1.5-2 \times 10^6$
100 mm dish	$3.5-4 \times 10^6$