



Thawing Hybridoma Cells

Day 1:

- 1) take cryovial from liquid nitrogen storage and quickly place on dry ice
- 2) pipet 1 ml of media pre-warmed to 37°C into a 10 cm Petri dish
- 3) immediately place cryovial in 37°C water bath and thaw watching closely (about 30-60 sec, just until cells are mostly thawed)
- 4) pipet the cells into the Petri dish containing 1 ml media
- 5) wait about 20-30 sec, then add 2 ml of media
- 6) wait about 20-30 sec, then add another 2 ml media
- 7) wait about 20-30 sec, then add another 2 ml media
- 8) add 8 ml media and place in 37°C (5% CO₂) incubator

Day 2:

Check cells under microscope to make sure they look healthy (cells should be round, bright, and clear)

If cells look healthy and are sufficiently confluent (50% confluent or greater), split 1:3 into new Petri dishes with media.

If most cells do not look healthy, harvest cells and centrifuge at 500 x g, aspirate media, resuspend in 15 ml of fresh pre-warmed media and replat into a fresh 10 cm Petri dish.

Day 3:

Make sure cells are healthy and split 1:3 into new Petri dishes with media

Media:

10% FetalClone I (ThermoFisher Cat# SH3008003)

1X PenStrep (ThermoFisher Cat# 15140122)

1X NEAA (ThermoFisher Cat# 11140050)

In DMEM (ThermoFisher Cat# 11965092)

Note that there are many suitable alternative protocols (e.g., Growing Hybridomas, Chapter 8. Antibodies: A Laboratory Manual. Second Edition. Greenfield, E. A. 2014. Cold Spring Harbor Laboratory Press.)