

# GENOTYPING PROTOCOL

## MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

[mmrrc@ucdavis.edu](mailto:mmrrc@ucdavis.edu)

530-754-MMRRC

**Protocol Name:** MMRRC 71890 C57BL/6N-Ace2em1(ACE2)Mbp Tmprss2em1(TMPRSS2)Mbp Plscrem1Mbp/Mmucd

**Protocol:** GoTaq® G2 Colorless Master Mix(Promega)

Reagent/Constituent	Volume (µL)
Water	4.0
GoTaq® G2 Colorless Master Mix,2X	7.5
Primer 1. (stock concentration is 20µM) IVF	0.5
Primer 2. (stock concentration is 20µM) IVR	0.5
Primer 3. (stock concentration is 20µM) kiR	0.5
DNA (example) extracted w/ "Qiagen DNeasy columns or other similar silica based kits"	1.5
<b>TOTAL VOLUME OF REACTION:</b>	
	<b>15.0 µL</b>

**Comments on protocol:**

- Protocol may work with other DNA extraction methods.

**Strategy:**

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting <span style="float: right;">HOT START? <input type="checkbox"/></span>	94	2:00	<b>1x</b>
2. Denaturation	94	0:10	
3. Annealing <span style="float: right;">steps 2-3-4 cycle in sequence</span>	65 (↓1°C/cycle)	0:30	<b>10x</b>
4. Elongation	68	2:00	
5. Denaturation	94	0:15	
6. Annealing <span style="float: right;">steps 5-6-7 cycle in sequence</span>	55	0:30	<b>25x</b>
7. Elongation	68	2:00 (↑20sec/cycle)	
8. Finish	4	∞	n/a

**Primers:**

**Electrophoresis Protocol:**

Name	Nucleotide Sequence (5' - 3')	Agarose: 1.5%	V: 90		
1. C19-Ace2-IVF	CTGTTTACATATCTGTCCTCTCCAGG	Estimated Running		90 min.	
2. C19-Ace2-IVR	GCTACAGAGGCAGTCACTCATCCTC	<b>Primer Combination</b>	<b>Band (bp)</b>	<b>Genotype</b>	
3. C19-hAce2-kiR	CCTCAGATCTCCAGCTTTCCCAA	1 & 2	522	wildtype	
		1 & 3, 1 & 2	723, 3757	mutant	



Name	Nucleotide Sequence (5' - 3')	Agarose: 1.5%	V: 90		
1. C19-Tmprss2-IVF	AGGTTCTCTGTACCTCAGAGGAGGA	Estimated Running		90 min.	
2. C19-Tmprss2-IVR	CCTGTCTCACCCCTTTCCAACATAACC	<b>Primer Combination</b>	<b>Band (bp)</b>	<b>Genotype</b>	
3. C19_Tmprss2-kiR	ACCTGAGGAGTCGCACTCTATCC	1 & 2	721	wildtype	
		1 & 3, 1 & 2	643, 4001	mutant	

Name	Nucleotide Sequence (5' - 3')	Agarose: 1.5%	V: 90		
1. CR-Plscr1-IVF	TTGAGCCACCATATGATTGCTG	Estimated Running		90 min.	
2. CR-Plscr1-IVR	GACTCTAGACAGCCACTGATAAGG	<b>Primer Combination</b>	<b>Band (bp)</b>	<b>Genotype</b>	
3. Post PCR digest w/BtsI-v2		1 & 2	437		
		Digest	159/278	wildtype	

# GENOTYPING PROTOCOL

## MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

[mmrrc@ucdavis.edu](mailto:mmrrc@ucdavis.edu)  
530-754-MMRRC




\*\*\*

	ACE2	TMPRSS2
#119	wt	wt
#141	hom	het
#148	hom	hom
#153	hemi	het

\*\*\* Note: Separate WT and Mut rxn required for Qiaxcel imaging (TMPRSS2).

Large mutant band is not observed with these PCR protocols.



Bts1-v2 Restriction Enzyme Digest (NEB)

Bts1-v2 Order R0667

### Restriction Enzyme Digestion


#### Steps

1. Set up reaction as follows:

COMPONENT	50 µl REACTION
DNA	1 µg
10X rCutSmart Buffer	5 µl (1X)
Bts1-v2	1.0 µl (10 units)†
Nuclease-free Water	to 50 µl

2. Incubate at 37°C for 5–15 minutes as Bts1-v2 is Time-Saver qualified.

Bts1-v2 Restriction Enzyme Digest (NEB)



	Plscr1_Q271Y
#249	wt
#250	hom
#251	het
#252	wt