

GENOTYPING BY PCR PROTOCOL

MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

mmrrc@ucdavis.edu

530-754-MMRRC

KK - 1May2019:

B6.129X1(SJL)-Scn1a^{tm2.1Wac}/Mmucd

MMRRC: 041829-UCD

Protocol:

Reagent/Constituent	Volume (µL)
Water	5.6
GoTaq® G2 Colorless Master Mix, 2X	7.5
Primer 1. (stock concentration is 20µM)	0.45
Primer 2. (stock concentration is 20µM)	0.45
DNA (example) extracted w/ "Qiagen DNeasy columns or other similar silica based kits"	1.0
TOTAL VOLUME	

Comments on protocol:

- ✓ Protocol may work with other DNA extraction methods.
- ✓ Use Touch-Down cycling protocol-first 10 cycles anneal at 65°C decreasing in temperature by 1.0°C; next 30 cycles anneal at 55°C.

Strategy:

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting HOT START? <input type="checkbox"/>	94	2:00	1x
2. Denaturation	94	0:10	
3. Annealing steps 2-3-4 cycle in sequence	65 (↓1°C/cycle)	0:30	
4. Elongation	68	2:00	
			Go to 2 (9x)
5. Denaturation	94	0:15	
6. Annealing steps 5-6-7 cycle in sequence	55	0:30	
7. Elongation	68	2:00 (+20s/cycle)	
			Go to 5 (24x)

Primers:

Name	Nucleotide Sequence (5' - 3')	Electrophoresis Protocol:		
1. 41829-F	CTTGATGTGTTGAAATTCAC	Argarose: 1.5% V: 90		
2. 41829-R	TATAGAGTGTTTAATCTCAAC	Estimated Running Time: 90 min.		
		Primer Combination	Band (bp)	Genotype
		1 & 2	1019	mutant
		1 & 2	846	wildtype

