

GENOTYPING BY PCR PROTOCOL FORM

MUTANT MOUSE REGIONAL RESOURCE CENTER: UC DAVIS

Investigator/PI: Bruce Beutler, M.D. Address: _____ Contact: Eva Moresco

Telephone: 858-784-1000 x3382 FAX: 858-784-8444 email: evamarie@scripps.edu

DNA Extraction Method: NaOH _____ Proteinase K X Other _____

Protocol: NAME OF PCR: MMRRC Strain #31704-UCD, C57BL/6J-Kcnj8^{sos} (SOS)

Reagent/ Constituent	Volume (uL)
DNA Sample	1
10x Buffer (contains 15mM MgCl2)	5
dNTPs (stock concentration is 25mM)	0.4
Primer 1 (stock concentration is 20 uM)	1
Primer 2 (stock concentration is 20 uM)	1
Taq Polymerase	2.5
water	39.1
TOTAL VOLUME OF REACTION:	50 uL

Comments on protocol (e.g., different concentration of MgCl₂, etc): Use primers 1 & 2 for amplification and 3 & 4 for sequencing.

Strategy:

Steps	Temp (°C)	Time (min)	# of Cycles
1. Initiation/Melting HOT START?...CHECK HERE []	94	2	1
2. Denaturation	94	30 seconds	
3. Annealing	56	30 seconds	30
4. Elongation		1 min	
5. Amplification (i.e., 72°C, 10 min)	72	7	1
6. Finish (i.e., 4°C, indefinite)	4	n/a	n/a

Primers:

Primer Name	Nucleotide Sequence (5' - 3')
1: Kcnj8_sos_F	GCACTCGCTGACCTGACATCATTC
2: Kcnj8_sos_R	TGACCTGACTACAGAGACTGCTGG
3: Kcnj8_sos_seqF	CCTGTTGATAATCCCATCGAGAG
4: Kcnj8_sos_seqR	TACAGAGACTGCTGGGAAGAAAAC

Electrophoresis Protocol:

% Agarose: _____ V : _____

Estimated Running Time (min): _____

Primer combination	Band (bp)	genotype
(i.e. 1&2)	1250	
(i.e. 3&4)		
(i.e. 1&2&3)		

Mutation site (red) and flanking sequence:

WT ccgaggagat**t**ccagtgggg
 sos ccgaggaga**c**ccagtgggg