

GENOTYPING PROTOCOL

MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

mmrrc@ucdavis.edu

530-754-MMRRC

Protocol Name: CR10538 Alg13 N107S (h.N107S) HDR Stock #: 66957

Reagent/Constituent	Volume (μL)
QuantiTect Multiplex PCR Master Mix Cat No./ID 204541	5.0
Water	3.4
Target Probe mix	0.3
-21 μM Mutant Forward Primer	
-21 μM Mutant Reverse Primer	
-7 μM Mutant probe	
TCRD (endogenous control) mix	0.3
-21 μM TCRD Forward primer	
-21 μM TCRD Reverse Primer	
-7 μM TCRD probe	
Sample	1.0
TOTAL VOLUME OF REACTION:	10.00 μL

Comments on protocol:

Protocol may work with other DNA extraction methods. Reference: ABI User Bulletin #2 and #5 (updated 10/2001) for multiplex in same tube and validation of each assay to match relative efficiencies of reference and target primer/probe combinations. Also reference: Rapid and accurate determination of zygosity in transgenic animals by real-time quantitative PCR. TransRes (2002).

Strategy:

Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting	95	15:00	1
2. Denaturation	95	0:30	40x
3. Annealing/Elongation	60	1:00	40x
4. To step 2 for 40 cycles			

Primers:


Name	Nucleotide Sequence (5' - 3')
1. TCRD Forward Primer	CAGACTGGTTATCTGCAAAGCAA
2. TCRD Reverse Primer	TCTATGCCAGTTCCAAAAACATC
3. TCRD MGB VIC Probe	VIC-ATTATAACGTGCTCCTGG-MGB
4. TM_Alg13-F	TCTGGAGAAAGGCAAACCAC
5. TM_Alg13-R	CTTTGTGCAACTGCTTTGCA
6. Alg13 BHQ-1 FAM Probe	Fam-AA-pdC-AG-pdC-pdC-A-pdU-pdC-AA-pdU-pdU-pdU-GAA-pdC-pdU-pdU-G-pdC-BHQ-1

Allele Description: The mouse N107S model was created using optimized CRISPR Cas9 KI technology utilizing Ribonucleoprotein (RNP) in the presence of a synthetic single strand DNA repair template harboring the desired KI. Zygotes were electroporated and subsequent progeny were screened for the presence of the correctly targeted allele via homology directed repair (HDR). The KI nucleotides are 11 and 12 bp from the cleavage site, and silent PAM mutations were engineered into the PAM site and ssODN to prevent cleavage of the KI allele by Cas9. Key progeny were sequence confirmed.

ctgtttggagagcttgagaaaggcaaaccttggttagttgtaaatgaaagttaatgaacaGCcatcaatttgaaCtTgcaaagcagttgcacaaagaaggccatctctttattgtac

WT AAT > AGC KI

Sample	ΔCt	Genotype
Alg13-ntc		No Rxn
Alg13-WT	16.65	WT
CR10538-25	1.25	Het
CR10538-26	1.22	Het



Het

aATcatcagtttgaaTtG
WT
aGCcatcagtttgaaCtT
KI

Note: Homozygous animals will have a deltaCt around -2.2 (1 Ct less than heterozygous animals).

GENOTYPING PROTOCOL

MUTANT MOUSE RESOURCE & RESEARCH CENTER: UC DAVIS

mmrrc@ucdavis.edu

530-754-MMRRC

Alternative Genotyping Protocol

Standard PCR and Sequencing

Protocol: *GoTaq® G2 Colorless Master Mix(Promega)*

Reagent/Constituent	Volume (μL)
Water	5.0
GoTaq® G2 Colorless Master Mix,2X	7.5
Primer 1. (stock concentration is 20μM) IVF	0.5
Primer 2. (stock concentration is 20μM) IVR	0.5
DNA (example) extracted w/ "Qiagen DNeasy columns or other similar silica based kits"	1.5
TOTAL VOLUME OF REACTION:	15.0 μL

Comments on protocol:

- Protocol may work with other DNA extraction methods.
- When crossing Alg13 and Glt28d2 the Taqman picks up the "other" mutation. Sequencing will be needed if mixed.

Strategy:


Steps	Temp (°C)	Time (m:ss)	# of Cycles
1. Initiation/Melting HOT START? ☐	94	2:00	1x
2. Denaturation	94	0:10	
3. Annealing steps 2-3-4 cycle in sequence	65 (↓1°C/cycle)	0:30	10x
4. Elongation	68	2:00	
5. Denaturation	94	0:15	
6. Annealing steps 5-6-7 cycle in sequence	55	0:30	25x
7. Elongation	68	2:00 (↑20sec/cycle)	
8. Finish	4	∞	n/a

Primers:


Name	Nucleotide Sequence (5' - 3')	Electrophoresis Protocol:		
1. CR Alg13-IVF	ATAGACATTACTCCTGGGTGGTGGTG	Agarose: 1.5%	:	90
2. CR Agl13-IVR	ACCAGTGAAAACCCGAGTAGTACTAGTGA	Estimated	90	min.
		Primer	Band (bp)	Seq Primer
		1 & 2	922	CR Alg13-IVF

Sequencing across Exon 3 to verify the HDR using standard PCR and PCR purification is used to determine the presence of HDR.

Note: It is necessary to sequence all WT-sized bands to verify HDR when using only standard PCR.



Sample	Genotype
CR10538-23	Het
CR10538-24	Het
CR10538-25	Het
CR10538-26	Het



aATcatcagtttgaaTtG
WT
aGcCcatcagtttgaaCtT
KI